

Nasal Fix Plaster

Securement device

Intended use:

Cureplast Nasal Fix Plaster is a specially designed securement plaster intended to hold nasal tubes and oxygen catheters securely in position at the nose. Its practical design helps stabilize the tube or catheter and reduces the risk of accidental movement or displacement during use. The plaster provides reliable fixation while conforming comfortably to the nasal area, making it suitable for use in hospitals, clinics, and other healthcare settings.

Intended Use:

Cureplast Nasal Fix Plaster is intended for the securement of oxygen catheters, nasogastric tubes, nasoduodenal tubes, and other compatible nasal tubing at the nose, helping to prevent accidental displacement. For single use only.



product range:

	Reference#	Size,mm
Nasal Fix Plaster Large	NF-101227	80 x 87
Nasal Fix Plaster Medium	NF-101303	70 x 76
Nasal Fix Plaster Small	NF-101909	28 x 37

Material specifications:

Contains latex? No
Contains DEHP? No

Certification:

CE-marking? Yes
Product classification Class I (EU)
Quality system certified to ISO 13485

Strelisation

Sterile? No
Shelf life 3 years

Ref#	Shelf box	
	Qty-pcs	Dimensions-mm
NF-101227	100	115 x 70 x 140
NF-101303	100	115 x 70 x 140
NF-101909	100	60 x 70 x 100

Storage and handling:

